


1




Open box and unpack items


Contents;

- 1x USB Microscope with USB lead
- 1x Installation CD with full manual
- 1x Quick start guide

Full Installation manual on CD for following languages



2




Insert installation CD

Make sure the USB microscope is plugged in to your PC or notebook, then insert the mini CD supplied in the box.

If you have a slot loading CD-Rom player or the mini CD does not fit your drive please visit our website and you can download the latest drivers.

Please note there are 2 software programs to install.

3



Install driver software

After inserting the mini CD you should see the install screen above. You must install both the microscope driver and the micro capture software after plugging your USB microscope into your PC or notebook.


Once installed reboot your PC and then plug your USB microscope into your USB port. Then click on desktop icon called micro capture.

Full user manual can be found on the CD supplied


+

More exciting products from veho

Multi-Media Webcam



Mini Video Camera - Internet Edition



For more information visit our website www.veho-world.com